

## Wrist Tag

The Wrist Tag provides protection to wander-prone residents (RoamAlert® system) or pediatric patients (HALO® system). Each tag is identified by its own unique code. The system generates an Exit alarm if the tag is brought near to an exit protected by a Door Controller. The Wrist Tag can be worn on the individual's wrist or ankle and is attached using VeriChip's tear-proof band.

In addition to the standard blue Wrist Tag (part no. AR3TA02-00W), VeriChip offers the green Wrist Tag with pulse Technology (part no. AR3TA01-00W). This tag sends a regular signal that is monitored by the system. The system tracks the location of all tags, enabling you to locate any monitored individual at any time.



| Wrist Tag Specifications |  |
|--------------------------|--|
| <b>Output Frequency</b>  | 433.92 MHz   |
| <b>Input Frequency</b>   | 307 kHz  |
| <b>Battery Life</b>      | AR3TA02-00W: three years; AR3TA01-00W: one year          |
| <b>Temperature</b>       | 14° to 122°F (-10° to 50°C)                              |
| <b>Humidity</b>          | 0-95% RH @ 70°F (21°C), non-condensing                   |
| <b>Dimensions</b>        | Approximately 1.1 x 1.0 x 0.5 in (2.79 x 2.54 x 1.27 cm) |
| <b>Weight</b>            | Approximately 0.3 oz. (8.5 g)                            |

### Features

- Unique electronic serial number to identify each tag
- High capacity battery
- Smooth surface permits easy sterilization and prevents harboring of bacteria
- Attractive, high-impact case
- Waterproof